

Product Description:

The repliQa HiFi Assembly Mix simplifies the construction of
recombinant DNA through the simultaneous and seamless
assembly of multiple DNA fragments possessing terminal
regions of sequence overlap in a single, isothermal reaction.
Similar in principle to the Gibson Assembly® Method, the high
efficiency repliQa HiFi Assembly Mix is ideal for a range of
genetic engineering applications including routine molecular
cloning, site-directed mutagenesis, assembly of large
constructs for synthetic biology applications, and the
construction of diverse sequence libraries for directed
evolution studies. The concentrated, two-component format
allows flexibility in design of assembly reactions and
compatibility with less concentrated DNA samples. The
repliQa Mix has been optimized for use with a total input
quantity of DNA fragments in the range of 0.03 to 0.5 pmols.
The assembly of up to six DNA fragments is recommended,
though the repliQa Mix has been successfully used for more
complex assemblies.

Product Specifications 95190-010 Rev 01

Product Information		
repliQa™ HiFi Assembly Mix		
Part Number	95190-010	
Number of Reactions	10 Reactions*	
Reaction Size	20 μL	
Storage Temperature	-25°C to -15°C	
Lot Number	021291	
Reference Number	071717, 060117, D5081717	
Expiration Date	06/30/2018	

^{*}Variable depending on individual usage of HiFi Assembly Mix

Component Part Numbers:

84614 repliQa HiFi Assembly Enzyme Mix, 0.02mL 84615 repliQa 10X Assembly Reaction Buffer, 0.01mL 84618 repliQa Dpnl Restriction Enzyme, 0.05mL

Product Specifications		
95190		
Assay	DNA Assembly Quantification Assay	
Specification	≥15%	

Quality Control Analysis and Specifications:

The repliQa HiFi Assembly Mix is functionally tested for assembly of three 1-kb PCR fragments and quantitation of the 3 kb product.

Enzyme components were tested prior to assembly and found free of contaminating endonucleases and exonucleases. Enzyme purity was >95% as determined by SDS-PAGE and negligible *E.coli* genomic DNA contamination was confirmed by qPCR.

Limitations of Use

Quantabio and Ultraplex are registered trademarks of QIAGEN Beverly, Inc. Quanta Biosciences, qScript, Geltrack, ToughMix, PerfeCta, and Fastmix are registered trademarks of Quanta BioSciences Inc. Extracta, AccuStart, AccuMelt, and Accuvue are trademarks of Quanta BioSciences Inc. Applied Biosystems, StepOne, StepOnePlus and ROX are trademarks of Thermo Fisher Scientific and or its subsidiaries. Please contact QIAGEN-Beverly for more information.

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. SDS sheets relevant to this product are available upon request.